



Recombinant Human Leukemia Inhibitory Factor for the Enhanced Expansion of Neural Stem Cells

Application Note

Workflow Summary

Applying neural stem cells (NSC) as a cell therapy has created new hope for the treatment of various neural degenerative diseases as well as the reversal of neural damage once thought to be irreparable. Currently, clinical trials are being performed to evaluate NSC for the treatment of stroke, Parkinson's disease, spinal cord injuries, and many other indications [1]. As new potential applications emerge and new clinical trials begin, the need for scalable and safe sources of the foundational cell culture materials becomes essential to the future of this technology. The ability to quickly and reliably expand neural NSCs in culture has therefore been essential to the exploration of these novel, cell-based therapeutic approaches. The goal of this application note is to provide step-by-step instructions for using [InVitria's recombinant hLIF](#) in the efficient, safe, and scalable expansion of NSC in culture.

Introduction

Applying neural stem cells (NSC) as a [cell therapy](#) has created new hope for the treatment of various neural degenerative diseases as well as the reversal of neural damage once thought to be irreparable. Currently, clinical trials are being performed to evaluate NSC for the treatment of stroke, Parkinson's disease, spinal cord injuries, and many other indications [1]. As new potential applications emerge and new clinical trials begin, the need for scalable and safe sources of the foundational cell culture materials becomes essential to the future of this technology. The ability to quickly and reliably expand neural NSCs in culture has therefore been essential to the exploration of these novel, cell-based therapeutic approaches.

The yield and reproducibility of NSC expansion is improved when appropriate signaling molecules are introduced to the culture medium [2]. For example, the cytokine known as Leukemia Inhibitory Factor (LIF) has been shown to activate the LIF receptor on the



surface of NSC, thereby enhancing the cells' capacity for self-renewal [3]. LIF is therefore considered a powerful tool for improving expansion of NSC *in vitro* [2,3]. Currently, recombinant human (rhLIF) for cell-culture applications is predominantly biomanufactured using *E. coli* as the host organism [4]. While this protein-expression system has helped advance NSC technology, concerns exist about the safety and unintended effects of contamination from trace levels of bacterial endotoxin. Furthermore, different methods of production are required to support the manufacturing scale required for commercial development of stem cell therapies like NSC.

InVitria has improved both the scale and quality of rhLIF by biomanufacturing the cytokine in a eukaryotic but completely animal-free and non-mammalian host [5]. A side-by-side comparison of InVitria's rhLIF and rhLIF isolated from *E. coli* for NSC expansion found that the two products produced equivalent decreases in doubling time. [6]. Furthermore, [InVitria's rhLIF product](#) was found to contain 500-fold lower bacterial endotoxin compared to *E. coli*-derived rhLIF [5].

The goal of this application note is to provide step-by-step instructions for using InVitria's recombinant hLIF in the efficient, safe, and scalable expansion of NSC in culture.



Materials

Cells

- ✓ H9-embryonic-derived Neural Stem Cells (FisherScientific Prod#: N7800100)

Growth Media Preparation

- ✓ DMEM/F-12 (Sigma-Aldrich Prod#: D6429)
- ✓ L-glutamine (Sigma-Aldrich Prod#: G7513)
- ✓ rhLIF Animal Component Free (InVitria#: 777LIF048)
- ✓ rhFGF-2 (Sigma-Aldrich Prod#: F0291)
- ✓ rhEGF (Sigma-Aldrich Prod#: E9644)
- ✓ N-2 Media Supplement (ThermoFisher Prod#: 17502001)

Plate Design

- ✓ Corning™ BioCoat™ Poly-L-Ornithine/Laminin Multiwell Plate (FisherScientific Prod#: 354659)
- ✓ DPBS (Sigma-Aldrich Prod#: D8537)
- ✓ EDTA (Sigma-Aldrich Prod#: 1233508)
- ✓ TrypLE Express Enzyme (FisherScientific Prod#: LS12604013)

Cell Counting

- ✓ Cell counter or hemocytometer

Equipment

- ✓ 37 °C Incubator at 5% CO₂



Protocol

Growth Media Preparation

1. To create the standard background media formulation for NSC expansion, supplement L-glutamine into DMEM/F12 for a final concentration of 2mM, add N-2 supplement for a concentration of 1%, and finally add FGF-2 & EGF to the media at a concentration of 10 ng/ml each. This formulation will then be supplemented with different sources of recombinant hLIF to produce the experimental conditions.
2. In separate media containers, supplement with recombinant hLIF. Add recombinant hLIF sourced from InVitria at a concentration 10 ng/ml.

Cell Plating & Expansion

1. In Poly-L-ornithine/laminin-coated plates, plate cells at a concentration of 0.2×10^5 viable cells/cm² in growth media described above.
2. Allow cells to grow in an incubator for 96-120 hours at 37 °C and 5% CO₂.
3. To harvest cells, add DPBS supplemented with TrypLE and 1 mM EDTA.

Cell Counting

1. Determine the concentration of cells/ml using a hemocytometer or cell counter in order to calculate doubling time.



Results and Discussion

Supplementing the NSC growth medium with recombinant hLIF significantly decreased the doubling times for NSC expansion *in vitro* (Fig. 1). Additionally, no significant difference in growth performance was observed between recombinant rhLIF from InVitria and rhLIF derived from *E. coli*. These results suggest that InVitria's non-mammalian expression system is capable of producing recombinant cytokines that are functionally comparable, yet more scalable and at higher quality required for clinical application.

With no change in performance, 500-fold lower levels of endotoxin, and unmatched scalability, InVitria's rhLIF provides clear benefits over bacterial expression rLIF in the culture and expansion of NSC.




Figure 1. H9-derived NSC were expanded for 96 hours. Data were collected from 3 separate experiments with 8 replicates from each experiment. The inclusion of InVitria's animal-component-free rhLIF and *E. coli*-derived rhLIF both produced significant decreases in doubling time of NSC compared to culture conditions without rhLIF supplementation ($P > 0.05^*$). There was no significant difference between the InVitria's rhLIF and *E. coli* derived rhLIF in terms of doubling-time performance ($P > 0.05$) [6].



References

1. Trounson A, Thakar RG, Lomax G, Gibbons D. Clinical trials for stem cell therapies. *BMC Med.* 2011;9:52. <https://doi.org/10.1186/1741-7015-9-52>
2. Buono KD, Vadlamuri D, Gan Q, Levison SW. Leukemia inhibitory factor is essential for subventricular zone neural stem cell and progenitor homeostasis as revealed by a novel flow cytometric analysis. *Dev Neurosci.* 2012;34(5):449-462. <https://doi.org/10.1159/000345155>
3. Pitman M, Emery B, Binder M, Wang S, Butzkueven H, Kilpatrick TJ. LIF receptor signaling modulates neural stem cell renewal. *Mol Cell Neurosci.* 2004;27(3):255-266. <https://doi.org/10.1016/j.mcn.2004.07.004>
4. Gearing DP, Nicola NA, Metcalf D et al. Production of Leukemia Inhibitory Factor in *Escherichia coli* by a Novel Procedure and Its Use in Maintaining Embryonic Stem Cells in Culture. *Bio/Technology.* 1989;7(11):1157-1161. <https://doi.org/10.1038/nbt1189-1157>
5. Youngblood BA, Alfano R, Pettit SC et al. Application of recombinant human leukemia inhibitory factor (LIF) produced in rice (*Oryza sativa* L.) for maintenance of mouse embryonic stem cells. *J Biotechnol.* 2014;172:67-72. <https://doi.org/10.1016/j.jbiotec.2013.12.012>
6. Alfano R, Youngblood BA, Zhang D, Huang N, MacDonald CC. Human leukemia inhibitory factor produced by the ExpressTec method from rice (*Oryza sativa* L.) is active in human neural stem cells and mouse induced pluripotent stem cells. *Bioengineered.* 2014;5(3):180-185. <https://doi.org/10.4161/bioe.28996>

Written by Randy Alfano, PhD, VP Product Development, InVitria

Dr. Alfano leads the product development team, where he uses his expertise in media design and optimization. In addition, he has specialized in the development and optimization of cell culture processes in mammalian cells for cancer biologics, in vivo animal models for metastatic prostate cancer, and was instrumental in developing thorough in vitro and in vivo models for immunostimulatory antigens. Dr. Alfano received his Ph.D. from Texas A&M Health Science Center College of Medicine.

For additional product or technical information, please contact us at

2718 Industrial Drive
Junction City, KS 66441

Phone: 1.800.916.8311
Email: Info@InVitria.com

