



2718 Industrial Drive  
Junction City, KS 66441  
Tel: 800-916-8311  
Fax: 785-238-5688  
[www.InVitria.com](http://www.InVitria.com)  
[info@InVitria.com](mailto:info@InVitria.com)

# ITSE Animal-Free

## Product Specifications

**Product Name:** ITSE Animal-Free  
**Product Number:** 777ITS032  
**Product Grade:** Cell culture  
**Product Form:** Liquid  
**Country of Origin:** USA

The product has been manufactured without the use of an animal host cell and without animal components. It is processed using equipment that is free of animal components. The packaging materials which contact the final product are also animal component free.

---

TEST	SPECIFICATION
Conductivity	13.7 – 16.7 mS/cm
pH	7.7 – 8.7
Sterility	Pass

---