



2718 Industrial Drive
Junction City, KS 66441
Tel 800-916-8311
Fax 785-238-5688
www.InVitria.com
info@InVitria.com

Product Specifications

Product Name: OptiPEAK T Lymphocyte Complete Media Kit
Kit Components: OptiPEAK T Lymphocyte (50X) Protein Supplement, OptiPEAK T Lymphocyte (10X) Media Supplement and OptiPEAK T Lymphocyte Base Media
Product Number: 777OPT069
Product Form: Liquid
Country of Origin: U.S.A.

The product has been manufactured without the use of an animal host cell and without animal components. It is processed using equipment that is free of animal components. The packaging materials which contact the final product are also animal component free.

TEST	SPECIFICATION
------	---------------

OptiPEAK T Lymphocyte (50X) Protein Supplement

pH	6.90 – 7.50
Osmolarity (mOsm/L)	300 – 330
Endotoxin	≤ 3.000 EU/ml
Sterility	Sterile

OptiPEAK T Lymphocyte (10X) Media Supplement

pH	6.50 – 7.20
Osmolarity (mOsm/L)	400 – 500
Endotoxin	≤ 1.000 EU/ml
Sterility	Sterile

OptiPEAK T Lymphocyte Base Media

pH	6.90 – 7.60
Osmolarity (mOsm/L)	260 – 290
Endotoxin	≤ 1.000 EU/ml
Sterility	Sterile

Kit Functionality in T Cells

Pass