



2718 Industrial Drive  
Junction City, KS 66441  
Tel 800-916-8311  
Fax 785-238-5688  
[www.InVitria.com](http://www.InVitria.com)  
[info@InVitria.com](mailto:info@InVitria.com)

## Product Specifications

**Product Name:** ITSE + A Animal-Free Solution  
**Product Number:** 777ITS092  
**Product Form:** Liquid  
**Country of Origin:** U.S.A.

The product has been manufactured without the use of an animal host cell and without animal components. It is processed using equipment that is free of animal components. The packaging materials which contact the final product are also animal component free.

---

TEST	SPECIFICATION
pH	7.00 – 8.50
Conductivity	5.0 – 20.0 mS/cm
Sterility	Sterile

---